



Anti-TMEM14E (aa 32-79) polyclonal antibody (DPABH-17801)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM14E (transmembrane protein 14E) is a protein-coding gene.
Immunogen	antigen sequence corresponding to amino acids 32-79 (KVGSVQSPSA GFFFSELAGL DASQPSRNPKEHLSSPVYIW DLARYYAN) of Human TMEM14E.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [TMEM14E transmembrane protein 15E \[Homo sapiens \]](#)

Official Symbol	TMEM14E
Synonyms	TMEM14E; transmembrane protein 14E;
Entrez Gene ID	645843
Protein Refseq	NP_001116700.1
UniProt ID	Q6UXP3