



Anti-STARD9 (aa 3020-3138) polyclonal antibody (DPABH-17811)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	antigen sequence, corresponding to internal sequence amino acids 3020-3138 of Human STARD9.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	STARD9 StAR-related lipid transfer (START) domain containing 10 [Homo sapiens]
Official Symbol	STARD9

Synonyms	STARD9; StAR-related lipid transfer (START) domain containing 9; stAR-related lipid transfer protein 9; START domain containing 9; START domain-containing protein 9;
Entrez Gene ID	57519
Protein Refseq	NP_065810.2
UniProt ID	Q9P2P6
Function	ATP binding; lipid binding; microtubule binding; microtubule motor activity