



Anti-SFT2D1 (aa 1-34) polyclonal antibody (DPABH-17871)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in fusion of retrograde transport vesicles derived from an endocytic compartment with the Golgi complex.
Immunogen	amino acids 1-34 (MEKLRRVLSG QDDEEQGLTA QVLDASSLSF NTRL) of Human STF2D1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [SFT2D1 SFT2 domain containing 2 \[Homo sapiens\]](#)

Official Symbol	SFT2D1
Synonyms	SFT2D1; SFT2 domain containing 1; pRGR1; C6orf83; vesicle transport protein SFT2A; SFT2 domain-containing protein 1;
Entrez Gene ID	113402
Protein Refseq	NP_660152.1
UniProt ID	Q8WV19