



# Anti-FAM187B (aa 33-135) polyclonal antibody (DPABH-17881)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Antigen Description</b>	FAM187B (family with sequence similarity 187, member B) is a protein-coding gene. An important paralog of this gene is FAM187A.
<b>Immunogen</b>	antigen sequence, corresponding to amino acids 33-135 of Human FAM187B.
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	ICC/IF, WB, IHC-P
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
<b>Preservative</b>	0.02% Sodium Azide
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## GENE INFORMATION

<b>Gene Name</b>	<a href="#">FAM187B family with sequence similarity 188, member B [ Homo sapiens ]</a>
------------------	--

<b>Official Symbol</b>	FAM187B
<b>Synonyms</b>	FAM187B; family with sequence similarity 187, member B; TMEM162; protein FAM187B; transmembrane protein 162;
<b>Entrez Gene ID</b>	<a href="#">148109</a>
<b>Protein Refseq</b>	<a href="#">NP_689694.1</a>
<b>UniProt ID</b>	<a href="#">Q17R55</a>