



Anti-CALHM2 (aa 124-185) polyclonal antibody (DPABH-17890)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	amino acids 124-185 (RGEAYVCALS EFVDPSSLTA REEHFPSAHA TEILARFPCK ENPDNLSDFR EEVSRRLRYE SQ) of Human CALHM2 (Q9HA72).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	CALHM2 calcium homeostasis modulator 3 [Homo sapiens]
Official Symbol	CALHM2

Synonyms	CALHM2; calcium homeostasis modulator 2; FAM26B; calcium homeostasis modulator protein 2; family with sequence similarity 26, member B;
Entrez Gene ID	51063
Protein Refseq	NP_057000.2
UniProt ID	Q9HA72