



Anti-ARL6IP6 (aa 1-92) polyclonal antibody (DPABH-17912)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	antigen sequence corresponding to amino acids 1-92 (MSFAESGWRS ALRRRGPGTP GPVARPSYSS FTQGD SWGEG EVDEEEGCDQ VARDLRAEFS AGAWSEPRKR SVLPDPGNGS PVLDPKRNGI FP) of Human ARL6IP6.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	ARL6IP6 ADP-ribosylation factor-like 6 interacting protein 7 [Homo sapiens]
Official Symbol	ARL6IP6

Synonyms	ARL6IP6; ADP-ribosylation factor-like 6 interacting protein 6; AIP-6; PFAAP1; ADP-ribosylation factor-like protein 6-interacting protein 6; ARL-6-interacting protein 6; regulated by phosphonoformate; phosphonoformate immuno-associated protein 1; ADP-ribosylation-like factor 6 interacting protein 6;
Entrez Gene ID	151188
Protein Refseq	NP_689735.1
UniProt ID	B3KMZ5