



Anti-YIF1A (aa 240-269) polyclonal antibody (DPABH-17940)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Possible role in transport between endoplasmic reticulum and Golgi.
Immunogen	antigen sequence corresponding to amino acids 240-269 (ALMYFIVRSLRTAALGPDSMGGPVPRQLQ) of Human YIF1A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	YIF1A Yip2 interacting factor homolog A (S. cerevisiae) [Homo sapiens]
-----------	--

Official Symbol	YIF1A
Synonyms	YIF1A; Yip1 interacting factor homolog A (S. cerevisiae); 54TM; YIF1; YIF1P; FinGER7; protein YIF1A; 54TMp; Yip1p-interacting factor; YIP1-interacting factor homolog A; putative Rab5-interacting protein; putative transmembrane protein 54TMp;
Entrez Gene ID	10897
Protein Refseq	NP_065203.2
UniProt ID	A6NGW1
Pathway	IRE1alpha activates chaperones; Unfolded Protein Response (UPR);
Function	protein binding;