



Anti-YIPF3 (aa 302-350) polyclonal antibody (DPABH-17949)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May play a role in hematopoiesis.
Immunogen	antigen sequence, corresponding to amino acids 302-350 (GILDTLEGPN IPPIQRVPRD IPAMLPAARL PTTVLNATAK AVAVTLQSH) of Human YIPF3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [YIPF3 Yip1 domain family, member 4 \[Homo sapiens \]](#)

Official Symbol	YIPF3
Synonyms	YIPF3; Yip1 domain family, member 3; KLIP1; FinGER3; C6orf109; dJ337H4.3; protein YIPF3; killer lineage protein 1; natural killer cell-specific antigen KLIP1;
Entrez Gene ID	25844
Protein Refseq	NP_056203.2
UniProt ID	Q9GZM5