



Anti-EQTN (aa 215-294) polyclonal antibody (DPABH-17977)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Involved in acrosome biogenesis.
Immunogen	antigen: SCESQYSVNP ELATMSYFHP SEGVSDTSFS KSAESSTFLG TTSSDMRRSG TRTSESKIMT DIISIGSDNE MHENDESVTR, corresponding to amino acids 215-294 of Human AFAF
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	EQTN equatorin, sperm acrosome associated [Homo sapiens]
Official Symbol	EQTN
Synonyms	EQTN; equatorin, sperm acrosome associated; AFAF; SPACA8; C9orf11; equatorin; equatorin; sperm acrosome associated 8; acrosome formation associated factor; acrosome formation-associated factor;
Entrez Gene ID	54586
Protein Refseq	NP_001155057.1
UniProt ID	Q9NQ60