



Anti-ZNF532 (aa 1094-1202) polyclonal antibody (DPABH-18019)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in transcriptional regulation.
Immunogen	antigen, corresponding to amino acids 1094-1202 of Human ZNF532.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [ZNF532 zinc finger protein 533 \[Homo sapiens \]](#)

Official Symbol	ZNF532
Synonyms	ZNF532; zinc finger protein 532;
Entrez Gene ID	55205
Protein Refseq	NP_060651.2
UniProt ID	A0A024R283
Function	DNA binding; metal ion binding;