



Anti-SFXN5 (aa 223-251) polyclonal antibody (DPABH-18030)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Transports citrate. Potential iron transporter.
Immunogen	antigen, corresponding to amino acids 223-251 (ELEEGIDVLDSG NLVGSSKIAARHALLE) of Human SFXN5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [SFXN5 sideroflexin 6 \[Homo sapiens \]](#)

Official Symbol	SFXN5
Synonyms	SFXN5; sideroflexin 5; BBG-TCC; sideroflexin-5;
Entrez Gene ID	94097
Protein Refseq	NP_653180.1
UniProt ID	Q8TD22
Function	cation transmembrane transporter activity; citrate transmembrane transporter activity;