



Anti-EQTN (aa 35-182) polyclonal antibody (DPABH-18041)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Involved in acrosome biogenesis.
Specificity	DPABH-18041 shows specific reactivity against Isoform 1 of Human AFAF.
Immunogen	antigen, corresponding to amino acids 35-182 of Human AFAF.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	EQTN equatorin, sperm acrosome associated [Homo sapiens]
Official Symbol	EQTN
Synonyms	EQTN; equatorin, sperm acrosome associated; AFAF; SPACA8; C9orf11; equatorin; equatorin; sperm acrosome associated 8; acrosome formation associated factor; acrosome formation-associated factor;
Entrez Gene ID	54586
Protein Refseq	NP_001155057.1
UniProt ID	Q9NQ60