



Anti-TBC1D15 (aa 38-165) polyclonal antibody (DPABH-18071)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Acts as a GTPase activating protein for RAB7A. Does not act on RAB4, RAB5 or RAB6.
Immunogen	antigen, corresponding to amino acids 38-165 of Human TBC1D15.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [TBC1D15 TBC1 domain family, member 16 \[Homo sapiens \]](#)

Official Symbol	TBC1D15
Synonyms	TBC1D15; TBC1 domain family, member 15; RAB7-GAP; TBC1 domain family member 15; GAP for RAB7; GTPase-activating protein RAB7; Tre2/Bub2/Cdc16 domain family member 15; Tre-2-budding uninhibited by benzimidazole-cell division cycle 16 domain, domain family member 15;
Entrez Gene ID	64786
Protein Refseq	NP_001139685.2
UniProt ID	Q8TC07
Function	Rab GTPase activator activity; protein binding;