



Anti-ZNF516 (aa 610-714) polyclonal antibody (DPABH-18151)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in transcriptional regulation.
Immunogen	antigen sequence, corresponding to amino acids 610-714 of Human ZNF516.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [ZNF516 zinc finger protein 517 \[Homo sapiens \]](#)

Official Symbol	ZNF516
Synonyms	ZNF516; zinc finger protein 516; HsT287;
Entrez Gene ID	9658
Protein Refseq	NP_055458.1
UniProt ID	Q2YDX2
Function	DNA binding; metal ion binding;