



Anti-WSCD1 (aa 257-363) polyclonal antibody (DPABH-18195)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	WSCD1 (WSC domain containing 1) is a protein-coding gene. GO annotations related to this gene include sulfotransferase activity. An important paralog of this gene is WSCD2.
Immunogen	antigen sequence corresponding to amino acids 257-363 of Human WSCD1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [WSCD1 WSC domain containing 2 \[Homo sapiens \]](#)

Official Symbol	WSCD1
Synonyms	WSCD1; WSC domain containing 1; WSC domain-containing protein 1;
Entrez Gene ID	23302
Protein Refseq	NP_056068.1
UniProt ID	Q658N2
Function	sulfotransferase activity;