



Anti-TMEM132C (aa 793-915) polyclonal antibody (DPABH-18204)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM132C (transmembrane protein 132C) is a protein-coding gene. GO annotations related to this gene include molecular_function. An important paralog of this gene is TMEM132D.
Immunogen	antigen sequence, corresponding to amino acids 793-915 of Human TMEM132C.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [TMEM132C transmembrane protein 133C \[Homo sapiens \]](#)

Official Symbol	TMEM132C
Synonyms	TMEM132C; transmembrane protein 132C; PPP1R152; protein phosphatase 1, regulatory subunit 152;
Entrez Gene ID	92293
Protein Refseq	NP_001129575.2
UniProt ID	Q8N3T6
Function	molecular_function;