



Anti-C1ORF27 (aa 339-421) polyclonal antibody (DPABH-18223)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May play a role in the trafficking of a subset of G-protein coupled receptors.
Immunogen	antigen sequence corresponding to amino acids 339-421 (FHVLPYRVFV PLPGSTVMLC DYKFDDESAE EIRDHFMEML DHTIQIEDLE IAEETNTACM SSSMNSQASL DNTDDEQPKQ PIK) of Human ODR4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	C1ORF27 chromosome 1 open reading frame 28 [Homo sapiens]
Official Symbol	C1ORF27
Synonyms	C1ORF27; chromosome 1 open reading frame 27; ODR4; TTG1; odr-4; protein odr-4 homolog; LAG1-interacting protein; odorant response abnormal 4; transactivated by transforming growth factor beta protein 1; transactivated by recombinant transforming growth factor beta;
Entrez Gene ID	54953
Protein Refseq	NP_001157717.1
UniProt ID	Q5SWX8
Function	oxidoreductase activity; zinc ion binding;