



Anti-PIGG (aa 175-315) polyclonal antibody (DPABH-18233)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Ethanolamine phosphate transferase involved in glycosylphosphatidylinositol-anchor biosynthesis. Transfers ethanolamine phosphate to the GPI second mannose.
Immunogen	antigen sequence, corresponding to amino acids 175-315 of Human PIGG.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	PIGG phosphatidylinositol glycan anchor biosynthesis, class G [Homo sapiens]
------------------	--

Official Symbol	PIGG
Synonyms	PIGG; phosphatidylinositol glycan anchor biosynthesis, class G; GPI7; LAS21; PRO4405; RLGS1930; GPI ethanolamine phosphate transferase 2; phosphatidylinositol-glycan biosynthesis class G protein;
Entrez Gene ID	54872
Protein Refseq	NP_001120650.1
UniProt ID	Q5H8A4
Pathway	Glycosylphosphatidylinositol(GPI)-anchor biosynthesis; Metabolism of proteins; Post-translational protein modification.
Function	CP2 mannose-ethanolamine phosphotransferase activity; phosphotransferase activity, for other substituted phosphate groups;