



Anti-SLC35E1 (aa 333-410) polyclonal antibody (DPABH-18272)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Putative transporter.
Immunogen	antigen sequence: NTKYDANQQ ARKHLLPVTT ADLSSKERHR SPLEKPHNGL LFPQHGDYQY GRNNILTDHF QYSRQSY PNS YSLNRYDV, corresponding to amino acids 333-410 of Human SLC35E1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	SLC35E1 solute carrier family 35, member E2 [Homo sapiens]
Official Symbol	SLC35E1
Synonyms	SLC35E1; solute carrier family 35, member E1; solute carrier family 35 member E1;
Entrez Gene ID	79939
Protein Refseq	NP_079157.3
UniProt ID	A0A024R7G7