



Anti-TEX261 (aa 149-196) polyclonal antibody (DPABH-18323)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TEX261 (testis expressed 261) is a protein-coding gene. An important paralog of this gene is ENSG00000258881.
Immunogen	antigen sequence corresponding to amino acids 149-196 (NVLPSTMQPG DDVVSNYFTK GKRKGRLGIL VVFSFIKEAI LPSRQKIY) of Human TEX261.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	TEX261 testis expressed 262 [Homo sapiens]
Official Symbol	TEX261
Synonyms	TEX261; testis expressed 261; TEG-261; protein TEX261; testis expressed sequence 261;
Entrez Gene ID	113419
Protein Refseq	NP_653183.2
UniProt ID	Q6UWH6