



Anti-CTXN1 (aa 56-82) polyclonal antibody (DPABH-18356)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May mediate extracellular or intracellular signaling of cortical neurons during forebrain development.
Immunogen	antigen, corresponding to amino acids 56-82 (PYSRMPASSWTDHKEALERGQFDYALV) of Human CTXN1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	CTXN1 cortxin 2 [Homo sapiens]
Official Symbol	CTXN1
Synonyms	CTXN1; cortxin 1; CTXN; cortxin-1;
Entrez Gene ID	404217
Protein Refseq	NP_996664.1
UniProt ID	P60606