



Anti-RAPH1 (aa 3-75) polyclonal antibody (DPABH-18367)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Mediator of localized membrane signals. Implicated in the regulation of lamellipodial dynamics. Negatively regulates cell adhesion.
Immunogen	antigen, corresponding to amino acids 3-75 (QLSDEEIDHG AEEDSDKEDQ DLDKMFGAWL GELDKLTQSL DSDKPMEPVK RSPLRQETNM ANFSYRFSIY NLN) of Human RAPH1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	RAPH1 Ras association (RalGDS/AF-6) and pleckstrin homology domains 2 [Homo sapiens]
Official Symbol	RAPH1
Synonyms	RAPH1; Ras association (RalGDS/AF-6) and pleckstrin homology domains 1; LPD; RMO1; PREL2; PREL-2; ALS2CR9; ALS2CR18; RalGDS/AF-6; ras-associated and pleckstrin homology domains-containing protein 1; lamellipodin; proline-rich EVH1 ligand 2; amyotrophic lateral sclerosis 2 (juvenile) chromosome region, candidate 9; amyotrophic lateral sclerosis 2 (juvenile) chromosome region, candidate 18; amyotrophic lateral sclerosis 2 chromosomal region candidate gene 9 protein; amyotrophic lateral sclerosis 2 chromosomal region candidate gene 18 protein;
Entrez Gene ID	65059
Protein Refseq	NP_976241.1
UniProt ID	Q70E73