



Anti-ZDHHC23 (aa 1-83) polyclonal antibody (DPABH-18378)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in NOS1 regulation and targeting to the synaptic membrane.
Immunogen	N terminal amino acids 1-83 (MTQKGSMKPV KKKKTEEP EPLCCCEYID RNGEKNHVAT CLCDCQDLDE GCDRWITCKS LQPETCERIM DTISDRLRIP WLR) of Human ZDHHC23.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [ZDHHC23 zinc finger, DHHC-type containing 24 \[Homo sapiens \]](#)

Official Symbol	ZDHHC23
Synonyms	ZDHHC23; zinc finger, DHHC-type containing 23; NIDD; palmitoyltransferase ZDHHC23; DHHC-23; probable palmitoyltransferase ZDHHC23; zinc finger, DHHC domain containing 23; zinc finger DHHC domain-containing protein 23;
Entrez Gene ID	254887
Protein Refseq	NP_775841.2
UniProt ID	Q8IYP9
Function	protein-cysteine S-palmitoyltransferase activity; zinc ion binding;