



Anti-TMEM19 (aa 281-310) polyclonal antibody (DPABH-18412)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes the human homolog of the yeast Imp3 protein. The protein localizes to the nucleoli and interacts with the U3 snoRNP complex. The protein contains an S4 domain.
Immunogen	antigen sequence: YTGLDESTGMVVNSPTNKARHIAGKPILDN, corresponding to amino acids 281-310 of Human TMEM19.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	TMEM19 transmembrane protein 20 [Homo sapiens]
Official Symbol	TMEM19
Synonyms	TMEM19; transmembrane protein 19;
Entrez Gene ID	55266
Protein Refseq	NP_060749.2
UniProt ID	A0A024RBA1