



Anti-CCDC107 (aa 135-241) polyclonal antibody (DPABH-18423)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CCDC107 (coiled-coil domain containing 107) is a protein-coding gene.
Immunogen	Recombinant fragment corresponding to amino acids 135-241 of Human CCDC107
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59.98% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Store at -20°C.

GENE INFORMATION

Gene Name	CCDC107 coiled-coil domain containing 108 [Homo sapiens]
Official Symbol	CCDC107

Synonyms	CCDC107; coiled-coil domain containing 107; PSEC0222; coiled-coil domain-containing protein 107;
Entrez Gene ID	203260
Protein Refseq	NP_001182129.1
UniProt ID	Q8WV48