



Anti-PLBD2 (aa 448-565) polyclonal antibody (DPABH-18487)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Putative phospholipase.
Immunogen	Recombinant Protein Epitope Signature Tag, corresponding to amino acids 448-565 of PLBD2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59.98% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [PLBD2 phospholipase B domain containing 3 \[Homo sapiens \]](#)

Official Symbol	PLBD2
Synonyms	PLBD2; phospholipase B domain containing 2; P76; putative phospholipase B-like 2; PLB homolog 2; 76 kDa protein; LAMA-like protein 2; lamina ancestor homolog 2; phospholipase B-like 2 32 kDa form; phospholipase B-like 2 45 kDa form; phospholipase B domain-containing protein 2; mannose-6-phosphate protein associated protein p76;
Entrez Gene ID	196463
Protein Refseq	NP_001153199.1
UniProt ID	Q8NHP8
Function	hydrolase activity;