



Anti-GINM1 (aa 37-145) polyclonal antibody (DPABH-18503)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	GINM1 (glycoprotein integral membrane 1) is a protein-coding gene.
Immunogen	Recombinant Protein Epitope Signature Tag, corresponding to amino acids 37-145 of Human C6orf72.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59.9% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [GINM1 glycoprotein integral membrane 2 \[Homo sapiens \]](#)

Official Symbol	GINM1
Synonyms	GINM1; glycoprotein integral membrane 1; C6orf72; dJ12G14.2; glycoprotein integral membrane protein 1;
Entrez Gene ID	116254
Protein Refseq	NP_620140.1
UniProt ID	Q9NU53