



Anti-SLFN5 (aa 70-210) polyclonal antibody (DPABH-18569)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May have a role in hematopoietic cell differentiation.
Immunogen	antigen sequence, corresponding to internal sequence amino acids 70-210 of Human SLFN5.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [SLFN5 schlafen family member 6 \[Homo sapiens \]](#)

Official Symbol	SLFN5
Synonyms	SLFN5; schlafen family member 5;
Entrez Gene ID	162394
Protein Refseq	NP_659412.3
UniProt ID	Q08AF3
Function	ATP binding;