



Anti-FRRS1 (aa 187-336) polyclonal antibody (DPABH-18580)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Ferric-chelate reductases reduce Fe(3+) to Fe(2+) before its transport from the endosome to the cytoplasm.
Immunogen	antigen sequence, corresponding to internal sequence amino acids 187-336 of Human FRRS1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [FRRS1 ferric-chelate reductase 2 \[Homo sapiens \]](#)

Official Symbol	FRRS1
Synonyms	FRRS1; ferric-chelate reductase 1; SDR2; SDFR2; SDR-2; stromal cell-derived receptor 2; stromal cell derived factor receptor 2 homolog;
Entrez Gene ID	391059
Protein Refseq	NP_001013682.2
UniProt ID	Q6ZNA5
Function	ferric-chelate reductase activity; metal ion binding;