



Anti-GPR179 (aa 109-253) polyclonal antibody (DPABH-18590)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Orphan receptor.
Immunogen	antigen sequence, corresponding to internal sequence amino acids 109-253 of Human GPR179.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [GPR179 G protein-coupled receptor 180 \[Homo sapiens \]](#)

Official Symbol	GPR179
Synonyms	GPR179; G protein-coupled receptor 179; CSNB1E; GPR158L; GPR158L1; probable G-protein coupled receptor 179; GPR158-like 1; probable G-protein coupled receptor 158-like 1;
Entrez Gene ID	440435
Protein Refseq	NP_001004334.3
UniProt ID	A0A087X0K8
Function	G-protein coupled receptor activity;