



Anti-BEND3 (aa 15-77) polyclonal antibody (DPABH-18612)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	BEND3 (BEN domain containing 3) is a protein-coding gene.
Immunogen	antigen sequence: REELDMDVMR PLINEQNFDG TSDEEHEQEL LPVQKHYQLD DQEGISFVQT LMHLLKGNIG TGL, corresponding to internal sequence amino acids 15-77 of Human BEND3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	BEND3 BEN domain containing 4 [Homo sapiens]
Official Symbol	BEND3
Synonyms	BEND3; BEN domain containing 3; KIAA1553; BEN domain-containing protein 3;
Entrez Gene ID	57673
Protein Refseq	NP_001073919.1
UniProt ID	Q5T5X7