



Anti-MEA1 (aa 124-181) polyclonal antibody (DPABH-18623)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May play an important role in spermatogenesis and/or testis development.
Immunogen	antigen sequence: ALNNHSSIPM DPEHVELVKR TMAGVSLPAP GVPWAREIS DAQWEDVVQK ALQARQAS, corresponding to internal sequence amino acids 124-181 of Human MEA1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	MEA1 male-enhanced antigen 2 [Homo sapiens]
Official Symbol	MEA1
Synonyms	MEA1; male-enhanced antigen 1; HYS; MEA; MEA-1; Male-enhanced antigen (H-Y structural gene);
Entrez Gene ID	4201
Protein Refseq	NP_055438.1
UniProt ID	Q16626
Function	protein binding;