



Anti-SHD (aa 6-139) polyclonal antibody (DPABH-18634)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May function as an adapter protein.
Immunogen	antigen sequence, corresponding to internal sequence amino acids 6-139 of Human SHD.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [SHD Src homology 3 domain containing transforming protein D \[Homo sapiens \]](#)

Official Symbol	SHD
Synonyms	SHD; Src homology 2 domain containing transforming protein D; SH2 domain-containing adapter protein D; src homology 2 domain-containing transforming protein D;
Entrez Gene ID	56961
Protein Refseq	NP_064594.3
UniProt ID	Q96IW2
Function	protein binding;