



Anti-KIAA1462 (aa 1053-1186) polyclonal antibody (DPABH-18644)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	KIAA1462 (KIAA1462) is a protein-coding gene. Diseases associated with KIAA1462 include morbid obesity , and obesity.
Immunogen	antigen sequence, corresponding to internal sequence amino acids 1053-1186 of Human KIAA1462.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	KIAA1462 KIAA1463 [Homo sapiens]
Official Symbol	KIAA1462
Synonyms	KIAA1462; JCAD; junctional protein associated with coronary artery disease;
Entrez Gene ID	57608
Protein Refseq	NP_065899.1
UniProt ID	Q9P266