



Anti-KIAA0430 (aa 260-401) polyclonal antibody (DPABH-18665)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	This gene encodes a putative peroxisomal protein that appears to be conserved across Euteleostomi. In humans, it may be autoantigenic.
Specificity	DPABH-18665 shows specific reactivity against isoforms 1, 2 and 4 of Human LIMKAIN B1.
Immunogen	internal sequence amino acids 260-401 of Human LIMKAIN B1.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	KIAA0430 KIAA0431 [Homo sapiens]
Official Symbol	KIAA0430
Synonyms	KIAA0430; LKAP; MARF1; PPP1R34; meiosis arrest female protein 1; limkain b1; limkain-b1; meiosis arrest female 1; protein phosphatase 1, regulatory subunit 34;
Entrez Gene ID	9665
Protein Refseq	NP_001171927.1
UniProt ID	Q9Y4F3
Function	RNA binding; molecular_function; nucleotide binding;