



Anti-CCDC186 (aa 487-618) polyclonal antibody (DPABH-18686)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	antigen, corresponding to internal sequence amino acids 487-618 of Human C10orf118.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	CCDC186 coiled-coil domain containing 187 [Homo sapiens]
Official Symbol	CCDC186

Synonyms	CCDC186; coiled-coil domain containing 186; C10orf118; uncharacterized protein C10orf118; CTCL tumor antigen L14-2; CTCL tumor antigen HD-CL-01/L14-2;
Entrez Gene ID	55088
Protein Refseq	NP_060487.2
UniProt ID	Q7Z3E2