



Anti-LYSD3 (aa 114-214) polyclonal antibody (DPABH-18696)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	LYSD3 (LysM, putative peptidoglycan-binding, domain containing 3) is a protein-coding gene. An important paralog of this gene is LYSD4.
Immunogen	corresponding to internal sequence amino acids 114-214 of Human LYSD3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [LYSD3 LysM, putative peptidoglycan-binding, domain containing 4 \[Homo sapiens \]](#)

Official Symbol	LYSMD3
Synonyms	LYSMD3; LysM, putative peptidoglycan-binding, domain containing 3; lysM and putative peptidoglycan-binding domain-containing protein 3;
Entrez Gene ID	116068
Protein Refseq	NP_001273741.1
UniProt ID	Q7Z3D4