



Anti-DRAM2 (aa 233-266) polyclonal antibody (DPABH-18718)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Induces apoptotic cell death when co-expressed with DRAM1.
Immunogen	antigen sequence: ISLRVEANLH GLTLYDTAPC PINNERTRLL SRDI, corresponding to C terminal amino acids 233-266 of Human DRAM2.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, ICC/IF
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [DRAM2 DNA-damage regulated autophagy modulator 3 \[Homo sapiens \]](#)

Official Symbol	DRAM2
Synonyms	DRAM2; DNA-damage regulated autophagy modulator 2; PRO180; TMEM77; WWFQ154; DNA damage-regulated autophagy modulator protein 2; RP5-1180E21.1; transmembrane protein 77; damage regulated autophagy modulator 2;
Entrez Gene ID	128338
Protein Refseq	NP_848549.3
UniProt ID	A0A024R0B7