



Anti-PCNXL3 (aa 192-303) polyclonal antibody (DPABH-18729)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	PCNXL3 (pecanex-like 3 (Drosophila)) is a protein-coding gene. Diseases associated with PCNXL3 include breast cancer. An important paralog of this gene is PCNX.
Immunogen	antigen, corresponding to internal sequence amino acids 192-303 of Human PCNXL3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [PCNXL3 pecanex-like 4 \(Drosophila\) \[Homo sapiens \]](#)

Official Symbol	PCNXL3
Synonyms	PCNXL3; pecanex-like 3 (Drosophila); pecanex-like protein 3;
Entrez Gene ID	399909
Protein Refseq	NP_115599.2
UniProt ID	Q9H6A9