



Anti-SPATA3 (aa 1-81) polyclonal antibody (DPABH-18822)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	SPATA3 (spermatogenesis associated 3) is a protein-coding gene. Diseases associated with SPATA3 include cryptorchidism. GO annotations related to this gene include molecular_function.
Immunogen	antigen, corresponding to amino acids 1-81 (MKKVKKKRSE ARRHRDSTSQ HASSNSTSQQ PSPESTPQQP SPESTPQHSS LETTSRQPAF QALPAPEIRR SSCCLLSPDA N) of Human SPATA3.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	SPATA3 spermatogenesis associated 4 [Homo sapiens]
Official Symbol	SPATA3
Synonyms	SPATA3; spermatogenesis associated 3; TSARG1; spermatogenesis-associated protein 3; testis spermatocyte apoptosis-related protein 1; testis and spermatogenesis cell-related protein 1; testis and spermatogenesis cell apoptosis related gene 1;
Entrez Gene ID	130560
Protein Refseq	NP_620712.2
UniProt ID	Q8NHX4
Function	molecular_function;