



Anti-RWDD2B (aa 252-319) polyclonal antibody (DPABH-18833)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	RWDD2B (RWD domain containing 2B) is a protein-coding gene. Diseases associated with RWDD2B include monosomy 21. GO annotations related to this gene include molecular_function. An important paralog of this gene is RWDD2A.
Immunogen	antigen, corresponding to amino acids 252-319 of Human RWDD2B.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	RWDD2B RWD domain containing 3B [Homo sapiens]
Official Symbol	RWDD2B
Synonyms	RWDD2B; RWD domain containing 2B; GL011; C21orf6; RWD domain-containing protein 2B;
Entrez Gene ID	10069
Protein Refseq	NP_058636.1
UniProt ID	P57060
Function	molecular_function;