



Anti-HSCB (aa 160-232) polyclonal antibody (DPABH-18852)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May act as a co-chaperone in iron-sulfur cluster assembly in mitochondria.
Immunogen	antigen, corresponding to amino acids 160-232 (QFLIEIMEIN EKLAEESEA AMKEIESIVK AKQKEFTDNV SSAFEQDDFE EAKEILTKMR YFSNIEEKIK LKK) of Human HSCB.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	HSCB HscB mitochondrial iron-sulfur cluster co-chaperone [Homo sapiens]
-----------	---

Official Symbol	HSCB
Synonyms	HSCB; HscB mitochondrial iron-sulfur cluster co-chaperone; JAC1; HSC20; DNAJC20; dJ366L4.2; iron-sulfur cluster co-chaperone protein HscB, mitochondrial; J-type co-chaperone HSC20; dnaJ homolog subfamily C member 20; DnaJ (Hsp40) homolog, subfamily C, member 20; HscB iron-sulfur cluster co-chaperone homolog;
Entrez Gene ID	150274
Protein Refseq	NP_741999.3
UniProt ID	Q8IWL3
Pathway	Metabolism of proteins;
Function	chaperone binding; metal ion binding; molecular_function; protein binding