



Anti-DSCR4 (aa 1-107) polyclonal antibody (DPABH-18863)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	DSCR4 (Down syndrome critical region gene 4) is a protein-coding gene. Diseases associated with DSCR4 include gestational choriocarcinoma , and down syndrome. GO annotations related to this gene include molecular_function.
Immunogen	amino acids 1-107 of Human DSCR4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	DSCR4 Down syndrome critical region gene 5 [Homo sapiens]
Official Symbol	DSCR4
Synonyms	DSCR4; Down syndrome critical region gene 4; DCRB; DSCRB; Down syndrome critical region protein 4; Down syndrome critical region protein B;
Entrez Gene ID	10281
Protein Refseq	NP_005858.1
UniProt ID	P56555
Function	molecular_function;