



Anti-SLC38A6 (aa 213-250) polyclonal antibody (DPABH-18881)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Probable sodium-dependent amino acid/proton antiporter.
Immunogen	KWSIPCPLTL NYVEKGFQIS NVTDDCKPKL FHFSKESA internal sequence amino acids 213-250 of Human SLC38A6.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [SLC38A6 solute carrier family 38, member 7 \[Homo sapiens\]](#)

Official Symbol	SLC38A6
Synonyms	SLC38A6; solute carrier family 38, member 6; NAT-1; probable sodium-coupled neutral amino acid transporter 6; N-system amino acid transporter 1; N system amino acid transporter NAT-1; Na(+)-coupled neutral amino acid transporter 6;
Entrez Gene ID	145389
Protein Refseq	NP_001166173.1
UniProt ID	Q8IZM9