



Anti-SERINC4 (aa 358-426) polyclonal antibody (DPABH-18887)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Incorporates a polar amino acid serine into membranes and facilitates the synthesis of two serine-derived lipids, phosphatidylserine and sphingolipids.
Immunogen	EASYLAEVFG PLWIVKVYSY EFQKPSLCFC CPETVEADKG QRGGAARPAD QETPPAPPVQ VQHLSYNYS internal sequence amino acids 358-426 of Human SERINC4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	SERINC4 serine incorporator 5 [Homo sapiens]
Official Symbol	SERINC4
Synonyms	SERINC4; serine incorporator 4;
Entrez Gene ID	619189
Protein Refseq	NP_001244960.1
UniProt ID	A6NH21