



Anti-ERGIC1 (aa 52-125) polyclonal antibody (DPABH-18897)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Possible role in transport between endoplasmic reticulum and Golgi.
Immunogen	antigen sequence: VVNELYVDDP DKDSGGKIDV SLNISLPNLH CELVGLDIQD EMGRHEVGHI DNSMKIPLNN GAGCRFEGQF SINK, corresponding to internal sequence amino acids 52-125 of Human ERGIC1
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	ERGIC1 endoplasmic reticulum-golgi intermediate compartment (ERGIC) 2 [Homo sapiens]
Official Symbol	ERGIC1
Synonyms	ERGIC1; endoplasmic reticulum-golgi intermediate compartment (ERGIC) 1; NET24; ERGIC32; ERGIC-32; endoplasmic reticulum-Golgi intermediate compartment protein 1; ER-Golgi intermediate compartment 32 kDa protein; endoplasmic reticulum-golgi intermediate compartment 32 kDa protein;
Entrez Gene ID	57222
Protein Refseq	NP_001026881.1
UniProt ID	Q969X5
Function	protein binding;