



Rabbit Anti-Human ATP8B2 Polyclonal Antibody (DPABH-18919)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant Protein, antigen sequence: LVLAYKDLDEEYYEEWAERRLQASLAQDSREDRLASIYEEVENNML
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Antigen affinity purified
Conjugate	Unconjugated
Applications	IHC
Format	Liquid
Size	100 µl
Buffer	40% glycerol and PBS (pH 7.2).
Preservative	0.02% Sodium Azide
Storage	Store at +4°C for short term storage. Long time storage is recommended at -20°C. Gently mix before use. Optimal concentrations and conditions for each application should be determined by the user.

BACKGROUND

Introduction The protein encoded by this gene belongs to the family of P-type cation transport ATPases, and

to the subfamily of aminophospholipid-transporting ATPases. The aminophospholipid translocases transport phosphatidylserine and phosphatidylethanolamine from one side of a bilayer to another. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]

Keywords

ATP8B2; ATPase, aminophospholipid transporter, class I, type 8B, member 2; ATPID; phospholipid-transporting ATPase ID; ATPase, class I, type 8B, member 2; probable phospholipid-transporting ATPase ID; P4-ATPase flippase complex alpha subunit ATP8B2;

GENE INFORMATION

Entrez Gene ID

[57198](#)

UniProt ID

[P98198](#)
