



Anti-RASAL2 (aa 860-969) polyclonal antibody (DPABH-18942)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Inhibitory regulator of the Ras-cyclic AMP pathway.
Immunogen	antigen, corresponding to amino acids 860-969 of Human RASAL2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [RASAL2 RAS protein activator like 3 \[Homo sapiens \]](#)

Official Symbol	RASAL2
Synonyms	RASAL2; RAS protein activator like 2; NGAP; ras GTPase-activating protein nGAP; RAS protein activator-like 2; Ras protein activator like 1;
Entrez Gene ID	9462
Protein Refseq	NP_004832.1
UniProt ID	A8K2P3
Pathway	Ras signaling pathway;
Function	Ras GTPase activator activity; protein binding;