



Anti-CDPF1 (aa 8-113) polyclonal antibody (DPABH-18962)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CDPF1 (cysteine-rich, DPF motif domain containing 1) is a protein-coding gene.
Immunogen	antigen, corresponding to amino acids 8-113 of Human C22orf40
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [CDPF1 cysteine-rich, DPF motif domain containing 2 \[Homo sapiens \]](#)

Official Symbol	CDPF1
Synonyms	CDPF1; cysteine-rich, DPF motif domain containing 1; C22orf40; cysteine-rich DPF motif domain-containing protein 1; UPF0595 protein C22orf40;
Entrez Gene ID	150383
Protein Refseq	NP_997210.3
UniProt ID	Q6NVV7