



Anti-EDARADD (aa 110-211) polyclonal antibody (DPABH-18970)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Adapter protein that interacts with EDAR DEATH domain and couples the receptor to EDA signaling pathway during morphogenesis of ectodermal organs. Mediates the activation of NF-kappa-B.
Immunogen	antigen sequence: PTISDLLNDQ DLLDVIRIKL DPCHPTVKNW RNFASKWGMS YDELCFLEQR PQSPTLEFLL RNSQRTVGQL MELCRLYHRA DVEKVLRRWV DEEWPKRERG DP, corresponding to amino acids 110-211 of Human EDARADD
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	EDARADD EDAR-associated death domain [Homo sapiens]
Official Symbol	EDARADD
Synonyms	EDARADD; EDAR-associated death domain; ED3; EDA3; ECTD11A; ECTD11B; ectodysplasin-A receptor-associated adapter protein; crinkled homolog; EDAR-associated death domain protein; ectodysplasia A receptor associated death domain; ectodysplasin A receptor associated adapter protein;
Entrez Gene ID	128178
Protein Refseq	NP_542776.1
UniProt ID	Q8WWZ3